

TOPIC

Perfect and cost-effective
for healthcare facilities



TOPIC product range

TOPIC cleaning and disinfection appliances



The clean solution

TOPIC 10

For the clear-of-the-floor wall mounting (with or without slop sink)



The easy way of installing a cleaning and disinfection appliance. Space-saving and ready for connection.

Dimensions (mm)
w 500, d 450 or 588, h (appliance) 1,630

A practical and compact work unit consisting of the TOPIC 10 bedpan washer with a slop sink specially designed to accompany it.

Dimensions (mm)
w 1000, d 450 or 588, h (appliance) 1,630

TOPIC 20

Floor mounted appliances (with or without slop sink)



Floor mounted cleaning and disinfection appliance with plinth. Just place in position and connect. The ideal way of utilizing existing connections.

Dimensions (mm)
w 500, d 450 or 588, h 1,730

A work unit which is versatile in its functions, easy-to-install and compact. Consisting of the TOPIC 20 floor-mounted bedpan washer with plinth and a slop sink specially designed to accompany it.

Dimensions (mm)
w 1000, d 450 or 588, h (appliance) 1,730

TOPIC 40

The practical cabinet model fitting under a work surface



The TOPIC 40 cleaning and disinfection appliance is the perfect complement to the existing cabinet model. The appliance has a large work and storage surface. Equipped as standard with a lockable storage cabinet for containers for chemicals. Upstand available for an additional charge.

Dimensions (mm)
Appliance
w 900, d 600, h 900

The TOPIC 40 cleaning and disinfection appliance with built-in hand wash basin and small storage surface is a suitable and small versatile combined work unit. Equipped as standard with a lockable storage cabinet for containers for chemicals. Upstand available for an additional charge.

Dimensions (mm)
Appliance
w 900, d 600, h 1,000
Hand wash basin
w 340, d 370, h 150

TOPIC 40 E

Are you renovating or upgrading your facility?



The TOPIC 40 E is the perfect choice for quickly replacing existing appliances installed under work surfaces. Configurable to fit in all available recesses.

Dimensions (mm)
w 700, d 588, h 760 or 860

SAN 14 B

TOPIC 40 with slop sink and hand wash basin



The compact combined care unit for care needs in the dirty utility room, complete with work and storage surface, slop sink and hand wash basin. Equipped as standard with a lockable storage cabinet for containers for chemicals. Splash screen (optional).

Dimensions (mm)
w 1,400, d 600, h 1,000
Slop sink
h 700
Sink
Ø 364 mm
Hand wash basin
w 340, d 370, h 150

SAN 19 B

TOPIC 40 with slop sink, hand wash basin and cabinet



Combined care unit with a large work surface and a 500 mm wide cabinet. With work and storage surface as well as slop sink. Equipped as standard with a lockable storage cabinet for containers for chemicals. Splash screen (optional).

Dimensions (mm)
w 1,900, d 600, h 900
Slop sink
h 700
Sink
Ø 364 mm
Hand wash basin
w 340, d 370, h 150
Cabinet space
w 500

SAN 20 BW

TOPIC 40 with slop sink, hand wash basin, wash basin and cabinet



Combined care unit with a large 600 mm wide cabinet and a slop sink. Large work and storage surface. Equipped as standard with a lockable storage cabinet for containers for chemicals. Splash screen (optional).

Dimensions (mm)
w 2,000, d 600, h 1,000
Slop sink
h 700
Sink
Ø 364 mm
Hand wash basin
w 340, d 370, h 150
Wash basin
w 500, d 400, h 300
Cabinet space
w 600

SAN 24 W

TOPIC 40 with slop sink, wash basin and cabinet 1,000 mm



A combined care unit with a large 1,000 mm wide cabinet – a very popular choice. Equipped as standard with a lockable storage cabinet for containers for chemicals. Splash screen (optional).

Dimensions (mm)
w 2,400, d 600, h 900
Slop sink
h 700
Sink
Ø 364 mm
Wash basin
w 500, d 400, h 300
Cabinet space
w 1,000

SAN 29 BW

TOPIC 40 with slop sink, hand wash basin, wash basin and cabinet 1,500 mm



Large combined unit with an ideal length. All cleaning tasks are optimally combined into a single unit. The basins are ergonomically arranged. The table dimensions make working on this generously proportioned work-surface easy. Equipped as standard with a lockable storage cabinet for containers for chemicals. Splash screen (optional).

Dimensions (mm)
w 2,900, d 600, h 900
Slop sink
h 700
Sink
Ø 364 mm
Hand wash basin
w 340, d 370, h 150
Wash basin
w 500, d 400, h 300
Cabinet space
w 1,500



TOPIC bedpan washers are supplied with the innovative MIKE electronic control system as standard.

Pictures contain special equipment

Accessories and range of fixtures

In addition to the MEIKO cleaning and disinfection appliances, we can also supply you with a comprehensive range of accessories for safe and ergonomic work and hygienic storage of care utensils. 6 examples of practical solutions are illustrated below. We would be happy to develop a range of suggested solutions specially for you. Just give us a call and we would be pleased to help you at a moment's notice.



Drying and storage racks*
for disinfected urine bottles, bedpans and receptacles, with storage grid or drip tray. Available in widths of 500 mm to 1,400 mm.



Wall shelving*
with height-adjustable shelves available in all shelf lengths (depth 350 mm).



Cupboards with doors**
Cupboard units with hinged doors and one intermediate shelf. Available in widths of 500 mm, 600 mm, 800 mm, 1,000 mm, 1,200 mm and 1,500 mm.



Cupboards without doors
Cupboards for open storage with one intermediate shelf. Available in widths of 500 mm, 1,000 mm and 1,200 mm.



Sink unit*
For sluice units, complete with cupboards and drawers underneath; can be designed in various widths and heights.



Slop sink with integrated rim flushing*
Can be operated independently of the bed-pan washer (optional). Splash screen (optional).

* Pictures contain special equipment
** Optional: lamp-shade and lamp